




### COMMUNICATION CENTER

AOA Wisconsin  
13931 Spring Street  
Sturtevant, WI 53177

Phone: 1-800-262-5221  
Fax: 262-886-6879  
International: 262-886-1050



## LINGUAL INDIRECT Rx

Dr. \_\_\_\_\_ Acct # \_\_\_\_\_

Address \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Patient \_\_\_\_\_

Tel # \_\_\_\_\_ Fax # \_\_\_\_\_

E-Mail \_\_\_\_\_

Shipped \_\_\_\_\_ Placement Date \_\_\_\_\_  
(PLACEMENT DATE SHOULD BE 1-2 DAYS BEFORE ACTUAL INSERTION DATE)

PLEASE WRITE SPECIAL INSTRUCTIONS

#### Full Treatment

Supply AOA Band w/Bracket

STb (.018 ONLY)  
 Upper  Lower

Ormco Gen 7  
 .018  .022  
 Upper  Lower

Doctor is Supplying Brackets

Return weldable brackets for:

7 6	6 7
7 6	6 7

Indicate Molar Bracket choices on diagram  
TB (Twin Bracket) TT (Terminal Tube) HC (Hinge Cap)

CLASS® Traditional Set up  CLASS® STb Straight Wire Set up  
(Requires Diagnostic Set-up for Bracket Positioning)

Reset all Teeth  
 Reset only indicated teeth

7 6 5 4 3 2 1	1 2 3 4 5 6 7
7 6 5 4 3 2 1	1 2 3 4 5 6 7

Anterior Overjet  
 Ideal  
 Maintain  
 Set to \_\_\_\_ mm

Anterior Root Torque  
Upper Lower  
 Maintain  Maintain  
 Increase  Increase  
 Decrease  Decrease

Arch Form  
Upper Lower  
 Maintain  Maintain  
 Constrict  Constrict  
 Widen  Widen

Anterior Overbite  
 Ideal  
 Increase \_\_\_\_ mm  
 Decrease \_\_\_\_ mm

Indicate interproximal reduction (stripping on the diagram)

TARG (Brackets Positioned directly onto Malocclusion Model Omit Set-up instructions for this service.)

Standard Values

Upper	Tip	Ht	Torque	Lower	Tip	Ht	Torque
Centrals	5°	-	14°	Centrals	0°	-	0°
Laterals	8°	-	7°	Laterals	0°	-	0°
Cuspids	10°	+5	-2°	Cuspids	10°	+5	-7°
Bicuspid	0°	-	-7°	1st Bicuspid	0°	-	-12°
Molars	0°	-	-10°	2nd Bicuspid	0°	-	-16°
				1st Molar	0°	-	-22°
				2nd Molar	0°	-	-27°

Wires Traditional Archform				Wires STb Straight Wire Archform			
	U	L		U	L		
Respond	.0155	<input type="checkbox"/>	<input type="checkbox"/> TMA	.016	<input type="checkbox"/>	<input type="checkbox"/>	Can only select if CLASS Straight Wire set-up is selected
	.0175	<input type="checkbox"/>		.0175x.0175	<input type="checkbox"/>	<input type="checkbox"/>	
				.017x.025	<input type="checkbox"/>	<input type="checkbox"/>	CuNiTi™
NiTi™	.016*	<input type="checkbox"/>	<input type="checkbox"/> St Steel	.014	<input type="checkbox"/>	<input type="checkbox"/>	.013
	.018*	<input type="checkbox"/>		.016	<input type="checkbox"/>	<input type="checkbox"/>	.016
CuNiTi™	.013++	<input type="checkbox"/>		.018	<input type="checkbox"/>	<input type="checkbox"/>	.016x.016
	.017x.017*	<input type="checkbox"/>		.016x.016	<input type="checkbox"/>	<input type="checkbox"/>	.018x.018
	.017x.025	<input type="checkbox"/>		.016x.022	<input type="checkbox"/>	<input type="checkbox"/>	TMA .0175x.0175
				.017x.025	<input type="checkbox"/>	<input type="checkbox"/>	
D-Rect .016x.022	<input type="checkbox"/>	<input type="checkbox"/>		.018x.025	<input type="checkbox"/>	<input type="checkbox"/>	St Steel .018x.018

++ Straight Length \*Mushroom Shaped


#### Limited Treatment Social 6

STb (.018 ONLY)  Doctor is Supplying Brackets

Upper  Include 1st bicuspids  
 Lower  Include 1st bicuspids

Upper Lower  .013 CuNiTi (included)  
   .016 NiTi  
   .016 TMA  
   .018 NiTi

Return Chairside Adhesive  Do Not Provide Adhesive



Circle teeth that will be banded. Place an "X" on teeth to be extracted or missing. Indicate overcorrections by using arrows.

**Laboratory Use Only**

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

IB \_\_\_\_\_

MAC \_\_\_\_\_

MAC \_\_\_\_\_

MAC \_\_\_\_\_

SPM \_\_\_\_\_

AAC \_\_\_\_\_

AAC \_\_\_\_\_

\_\_\_\_\_

SHIP \_\_\_\_\_

RECEIVING

OPEN \_\_\_\_\_

ID \_\_\_\_\_

DATE \_\_\_\_\_

ENTER \_\_\_\_\_

PULL \_\_\_\_\_

SHIPPING

PULL \_\_\_\_\_

LAYOUT \_\_\_\_\_

PACK \_\_\_\_\_

SHIP \_\_\_\_\_

CHECK \_\_\_\_\_