

LAB USE ONLY

**Splint Rx**

Dr. \_\_\_\_\_ Acct # \_\_\_\_\_

Address \_\_\_\_\_

City, State, Zip \_\_\_\_\_

Patient \_\_\_\_\_

Tel # \_\_\_\_\_ Fax # \_\_\_\_\_

E-Mail \_\_\_\_\_

Shipped \_\_\_\_\_ Placement Date \_\_\_\_\_

(PLACEMENT DATE SHOULD BE 1-2 DAYS BEFORE ACTUAL INSERTION DATE)

**Please Ship Extra**

- Pre-Paid Bags     Shipping Boxes     Prescription Sheets

**Submitted scans digitally through:**

- Lythos     iTero     3M     Box     Other \_\_\_\_\_

**(Submit digital Rx to: [digital.services@aolab.com](mailto:digital.services@aolab.com))**


**Special Instructions:**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**Splint Instructional Drawings**



**Bite Registration Instructions**

- Checked and approved by Dr. - Use as is     Checked by Dr. - Call if questions

**For digital scans, please provide bite parameters:**

**Open Vertical** \_\_\_\_\_ mm    **Advance** \_\_\_\_\_ mm

**Allesee Orthodontic Appliances - 13931 Spring Street - Sturtevant, WI 53177**  
**Phone 1-800-262-5221 - Fax 262-886-6879 - International 262-886-1050**

**Flat Plane Occlusal Splint**

- Upper     Lower  
 Full coverage with even contact of opposing arch. Counter model suggested with wax bite.

**Superior Repositioning Splint**

- Upper     Lower (Tanner splint)  
 Full coverage on upper arch with anterior ramp designed to provide protrusive and cuspid protection. Both upper and lower models with *centric relation bite required*. Bite should provide 1-2 mm minimum anterior opening for acrylic durability.

Mounted models for adjustable articulators are recommended.

**Anterior Repositioning Splint**

- Upper     Lower  
 Full coverage upper appliance with anterior ramp designed to "capture" the lower anterior teeth and hold the mandible in a forward/ anterior position. Lower appliance consist of posterior coverage indexed with upper lingual cusps connected with an anterior lingual band of acrylic. Forward upper and lower models. *Bite reflecting anterior position.*

**Roth Splint Instructions**

- Repositioning Splint     Anterior Repositioning Splint  
 2 clasps (standard)    Length of splint ramp \_\_\_\_\_ mm  
 4 clasps    Clasp type     Adams     Ball

**Materials**

- Splint Grade Acrylic  
 Pressure formed Biocryl base -  Add Acrylic     No Acrylic  
 Soft Vinyl appliance -  2mm     3mm  
 Dual Hardness - Soft inner surface with hard outer shell  
 Color     Clear    Color \_\_\_\_\_

**Bleaching Tray**

- Upper     Lower     Vinyl no Bleach Reservoirs     Vinyl with Bleach Reservoirs  
*Indicate Reservoirs on Diagram*

**Clasps - Not available with Soft Vinyl or Dual Hardness Materials**

- Mark location on diagram  
 Ball     Adams     Arrow     "C"

**Retention Splint**

Requires upper and lower models and construction wax bite.

Please Check Bite Registration Instructions

**Damon**

- Hard Pressure formed  
 Dual Hardness with Soft Liner  
 Elastic Silicone  
 **Starnes Bite Orthotic**

- MORA**     **GELB**

(deep occlusal indents)    (point contact)  
 Lower splint with posterior coverage connected with a heavy metal bar lingual of the incisors. Counter model is recommended with a wax bite reflecting desired occlusal acrylic thickness.

**Anterior Open Bite Splint - AOB**

Upper posterior coverage with transpalatal bars. Buccal hooks are incorporated. Counter model advised. This appliance is used with implant screws (not available though AOA), which provides accelerated intrusion of upper posterior quadrants.

**Anterior Open Bite Retainer**

Doctor Signature \_\_\_\_\_

