

# POSITIONER

The **Positioner** has been considered by many clinicians to be the most complete means of finishing orthodontic treatment. Combining modern materials and fabrication techniques with the experience of AOA's technical staff has advanced the positioner to a multipurpose appliance. Its reduced wearing time and light weight mean improved patient cooperation and a refined finish. Even if the patient doesn't wear the positioner as prescribed, the super elastic material will engage and correct alignment.



Our set-up technicians can process your gnathological cases on virtually any articulator system, including the High Post Panadent and the Sam III instruments.

Our positioner department has the skill to fabricate accurate and cost efficient "bench" set-ups and positioners. We also specialize in partial and diagnostic set-ups.

The majority of the positioner cases forwarded to AOA have the brackets reflected on the models. We remove the brackets and fabricate the set-up and positioner. This method allows the "immediate" placement of the positioner upon clinical bracket removal when teeth are most susceptible to guidance.

All of our positioners and mouthguards are custom trimmed to your exact specifications. Their thickness and height can be tailored to produce the perfect blend of patient comfort and appliance efficacy. Custom mouthguard information is obtainable by calling our Centralized Communication Center.

**AOA**

COMMUNICATION CENTER  
**800.262.5221**  
INTERNATIONAL 262.886.1050  
FAX 262.886.6879

Please visit us at  
[www.aoaaccess.com](http://www.aoaaccess.com)

